

Group A

Exp # 13: Fluorescence of Anthracene

2026/04/1

By Sigi Li & Yuchen Shi

1. Preparation of Solution

Given cyclohexene, iodoethane (98%) and anthracene-cyclohexane (1:1) solution.

* From calculations in advance

- ① 0.050M quencher solution: $100\mu\text{L} = 0.1\text{mL}$
- ② 0.100M quencher solution: $200\mu\text{L} = 0.2\text{mL}$
- ③ 0.150M quencher solution: $300\mu\text{L} = 0.3\text{mL}$
- ④ 0.200M quencher solution: $400\mu\text{L} = 0.4\text{mL}$

2. Measurement

(1) Emission Spectrum at $\lambda_{\text{ex}} = 340\text{nm}$

(2) Emission Spectrum at $\lambda_{\text{ex}} = 325\text{nm}$

$\Rightarrow$ Overlay

(4) Excitation Spectrum {

0 M quencher
0.050M quencher
0.100M quencher
0.150M quencher
0.200M quencher

$\Rightarrow$ Overlay

~~10/4~~